PULPDENT®


Flows-Rite™

Multi-Purpose Flowable Composite

 4 x 1.5gm Flows-Rite

R_x only

 **PULPDENT Corporation** 80 Oakland Street • Watertown, MA 02472 • U.S.A.
www.pulpdent.com • www.activabioactive.com • pulpdent@pulpdent.com

PULPDENT®

Flows-Rite™

PULPDENT®

Flows-Rite™

Multi-Purpose Flowable Composite

CE

0459

PRODUCT DESCRIPTION

Flows-Rite is a 68% filled, light cured flowable composite, veneer cement, liner, base and fissure sealant with the handling characteristics and esthetic properties most preferred by dentists. It is radiopaque, fluoride releasing and has a nominal particle size of 0.7 microns.

Flows-Rite is available in the following shades: A1, A2, A3, A3.5, B2, C2, and ShadeFusion. ShadeFusion is a unique chameleon shade. It blends shades and reflects adjacent hues to create a natural esthetic, and it is ideal for veneer cementation.

INSTRUCTIONS FOR USE

Bonding to Dentin

1. Clean, prepare and isolate tooth from salivary contamination. Place appropriate base/liner in deep cavities, if desired.
2. Apply Etch-Rite or Etch Royale phosphoric acid gel to dentin for 15 seconds. Rinse completely. Keep surface clean. Avoid saliva contamination.
3. For highest bond strength, apply DenTASTIC UNO or your preferred adhesive following manufacturer's instructions for use.
4. Place pink 20 gauge applicator tip on Flows-Rite syringe and using gentle pressure carefully start flow of composite. Do not force the syringe. Use a new applicator tip for each patient.
5. Apply material to desired area.
6. Light cure. Cures with all lights. Curing time for a halogen light with a minimum of 300 mW/cm² is 30 seconds. More powerful lights will cure faster.
7. Place cap on syringe immediately after use.

Bonding to Enamel

1. Clean enamel surface with oil-free cleaning paste. Rinse well and dry. Isolate treatment area from salivary contamination.
2. Apply Etch-Rite or Etch Royale phosphoric acid gel to enamel for 20 seconds. Rinse completely. Do not disturb this surface.
3. Dry etched surface with clean, uncontaminated air. Do not blow saliva onto etched surface. If etched surface does not appear chalky white, etch again for 20 seconds, rinse and dry. If salivary contamination occurs, etch again for 10 seconds, rinse and dry.

4. If using an adhesive bonding agent, apply adhesive following manufacturer's instructions. A slightly moist tooth surface may be indicated.
5. Place pink 20 gauge applicator tip on Flows-Rite syringe and using gentle pressure carefully start flow of composite. Do not force the syringe. Use a new applicator tip for each patient.
6. Apply material to desired area.
7. Light cure. Cures with all lights. Curing time for a halogen light with a minimum of 300 mW/cm² is 30 seconds. More powerful lights will cure faster.
8. Place cap on syringe immediately after use.

Note: Multi-dose syringes should either be encased in a fresh protective barrier for each patient or cleaned and disinfected between patients, as appropriate.

Caution

Uncured material may cause eye or skin irritation on contact. Dental professionals should wear safety glasses and surgical gloves.

Storage and Handling

- Store tightly sealed in original container at cool room temperature. Avoid direct light, extremes of temperature, contamination and sources of ignition.
- Shelf life of unopened product: 2 years from date of manufacture.
- Re-cap syringe immediately after use

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