# Spee-Dee™ Build-Up

## Multi-Purpose Core and Tooth Build-Up Resin Dual Cure • Fluoride Releasing • Hydrophilic Radiopaque

### PRODUCT DESCRIPTION

Spee-Dee Build-Up is a multi-purpose, dual cure, resin-based dental restorative and core build-up material. It is moisture friendly and has self-adhesive and self-etching properties that are compatible with all bonding agents. The material is fluoride releasing and radiopaque. It contains no Bisphenol A, no Bis-GMA, and no BPA derivatives.

Spee-Dee Build-Up is dense like dentin and therefore cuts like dentin. Rotary instruments transition smoothly from tooth to Spee-Dee, cut accurately, and do not ditch, gouge, gum up or chatter.

Spee-Dee Build-Up is used for core build-up, post cementation, and as a multi-purpose repair material for fractured cusps and broken down teeth, defects in crown preparations, lost partial or whole fillings, and chipped or broken crowns or bridges.

Posts and cores can be placed in one step ensuring a homogeneous, one-piece internal structure.

In restorations with retention form, Spee-Dee Build-Up's self-adhesive and self-etching properties may be sufficient for post and core procedures without the use of a bonding agent. However, to achieve highest bond strengths, etching and bonding agents are recommended.

### INDICATIONS

Multi-purpose, resin-based dental restorative and core build-up material

### CONTRAINDICATIONS

Not indicated for direct placement on the exposed pulp.

### HOW TO USE THE DOUBLE BARREL CARTRIDGE

1. Remove cap. If necessary bleed the cartridge so that base and catalyst are at the orifice of the barrels. Place a mixing tip on the double barrel cartridge.
2. To ensure an even mix of base and catalyst, dispense 2-3 mm of material onto a pad and discard this material.
3. Dispense evenly mixed material directly onto the tooth or into the restoration.
4. Discard mixing tip. Recap cartridge. Do not cross-contaminate base and catalyst.

### DIRECTIONS FOR USE

1. Remove the old restoration and/or decay, or prepare post space using post system of choice
2. To achieve highest bond strengths, etch and place bonding agent according to the manufacturer's instructions. Light cure. A dual cure bonding agent is recommended with self-cure and dual cure resin core materials.
3. Remove cap from double barrel cartridge. If necessary bleed the cartridge so that base and catalyst are at the orifice of the barrels.
4. Place a mixing tip on the double barrel cartridge and insert into an appropriate 1:1 cartridge dispenser. To ensure an even mix of base and catalyst, dispense 2-3 mm of material onto a pad and discard this material.
5. Dispense material into the preparation site, or dispense into the post hole and place the post. Light cure for 20 seconds. Build-up the core using "Crown Up" technique by placing material, curing, adding more material, curing, and repeating this until the desired build-up is accomplished. The material will also self-cure. Anaerobic self-cure setting time at 37°C is 3:15 minutes. Carefully recap cartridge. Do not cross-contaminate base and catalyst.
6. Prepare build-up to receive the restoration.

### WARNING

- For professional use by dental personnel only.
- Pulpdent Spee-Dee Build-Up may cause irritation to eyes and skin on contact. Take all necessary precautions to avoid contact with eyes and soft tissue. If materials contact eyes, rinse with copious amount of water. If materials contact skin, wash skin with soap and water. Seek medical attention if irritation persists. See MSDS for further information.
- Wear protective eyewear and gloves. Wash hands after use. Emergency eye wash fountain should be available.

### STORAGE

Store tightly capped in original container at cool room temperature. Avoid direct sunlight, extremes of temperature and cross-contamination between parts.

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**PULPDENT Corporation**
80 Oakland Street • Watertown, MA 02472 • U.S.A.