

PULPDENT CORPORATION

Press Release

For further information contact:

Pulpdent Corporation

80 Oakland Street • Watertown, MA • 02471-0780 USA

(617) 926-6666 / (800) 343-4342 / pulpdent@pulpdent.com

PULPDENT INTRODUCES SPEE-DEE™ BUILD-UP

WATERTOWN, MA: March 3, 2009— Pulpdent Corporation has released Spee-Dee Build-Up, a dual cure multi-purpose core and tooth build-up resin. Spee-Dee provides one-step post and core build-ups to ensure a homogenous, one-piece internal structure. It is also used for post cementation, vital tooth build-up for crown preps, and patching broken cusps, porcelain, lost fillings and other defects.

Pulpdent Spee-Dee Build-up simulates tooth structure and really does cut like dentin. Rotary instruments transition smoothly from tooth structure to Spee-Dee Build-Up and cut accurately without gouging, gumming up or chattering. Spee-Dee is self-adhesive and self-etching to dentin, making the use of etching and bonding agents optional.

Spee-Dee Build-Up represents another application of Pulpdent's unique moisture-tolerant resin technology. Like the revolutionary Embrace™ WetBond™ products, Spee-Dee behaves favorably in the moist oral environment. It is dual-cure, fluoride-releasing, and radiopaque.

For more information call 800-343-4342 or visit www.pulpdent.com.

###