

PULPDENT CORPORATION

Press Release

For further information contact:

Pulpdent Corporation

80 Oakland Street • Watertown, MA • 02471-0780 USA

(617) 926-6666 / (800) 343-4342 / pulpdent@pulpdent.com

Pulpdent Introduces Prep-Rite RC Canal Lubricant

WATERTOWN, MA: September 24, 2007 - Pulpdent Corporation has introduced Prep-Rite RC, a 17% viscous EDTA gel that is easily picked up on an endodontic file to facilitate preparation of the root canal. Prep-Rite is the latest addition to Pulpdent's range of EDTA products for endodontics.

Prep-Rite RC is a canal lubricant and chelating agent that has been buffered to a neutral pH. It effectively decalcifies canals and facilitates removing dentin from the canal walls, and it acts as a lubricant to minimize the binding and breaking of files. Prep-Rite is water-based and easily rinses out of the canal using traditional irrigating solutions.

Pulpdent also manufactures 17% EDTA liquid solution and File-Rite EDTA semi-gel, both of which are syringed directly into the canal. Prep-Rite RC is picked up on a file using the traditional technique.

Prep-Rite RC is available in a kit containing 4 x 5 gram syringes. For more information, call 800-343-4342 or visit www.pulpdent.com.

###