

PULPDENT CORPORATION

Press Release

For further information contact:

Pulpdent Corporation

80 Oakland Street • Watertown, MA • 02471-0780 USA

(617) 926-6666 / (800) 343-4342 / pulpdent@pulpdent.com

Flows-Rite Features ShadeFusion Chameleon Shade

WATERTOWN, MA: December 4, 2007— Flows-Rite Multi-Purpose Flowable Composite and Veneer Cement from Pulpdent features a unique chameleon shade, ShadeFusion, that blends shades and reflects adjacent hues to create a natural esthetic. Flows-Rite is a 68% filled, light-cured, flowable restorative, veneer cement, liner, base, and fissure sealant with the handling characteristics and esthetic properties most preferred by dentists.

Pulpdent Flows-Rite is radiopaque, fluoride releasing and has a nominal particle size of 0.7 microns. It cures with any light.

Flows-Rite is available in shades A1, A2, A3, A3.5, B2, B3, C2, C3, bleach white, and ShadeFusion. The chameleon ShadeFusion is ideal for incisal repairs and veneer cementation.

Flows-Rite is sold in packages containing 4 x 1.5 gm syringes of Flows-Rite plus 20 applicator tips. Pulpdent also offers an A-Shade Assortment Kit with one syringe each A1, A2, A3, and A3.5; and a Rainbow Assortment Kit with one syringe each A2, B2, C2, and ShadeFusion. For more information, call 800-343-4342 or visit www.pulpdent.com.

###