

PULPDENT CORPORATION

Press Release

For further information contact:

Pulpdent Corporation

80 Oakland Street • Watertown, MA • 02471-0780 USA

(617) 926-6666 / (800) 343-4342 / pulpdent@pulpdent.com

Pulpdent introduces Forendo Paste

WATERTOWN, MA: September 18, 2008— Pulpdent Corporation has introduced Forendo™ Paste, a non-setting, premixed calcium hydroxide paste with iodoform for root canal therapy. Forendo Paste is a temporary endodontic filling material that promotes healing and inhibits bacterial growth in the root canal.

Forendo Paste contains calcium hydroxide and iodoform in a silicone oil base. It has a pH greater than 12 and is radiopaque. The paste is bactericidal and bacteriostatic. It is recommended for routine use as a disinfecting root canal dressing between office visits to prevent flare-ups and for treatment of infected root canals, abscesses, periapical lesions, perforations, exudation, root resorption, traumatic injuries, and for apexification and apexogenesis. Forendo Paste can be used alone as a temporary canal disinfectant, in conjunction with gutta percha for permanent filling, and for pulpectomies and vital pulpotomies in deciduous teeth.

Forendo Paste is a new addition to Pulpdent's clinically proven range of calcium hydroxide pastes for root canal therapy that also includes Pulpdent Paste, TempCanal and Multi-Cal.

###