


Syringe Applicator Tips and Dispensers

The higher the gauge, the smaller the needle diameter.
The 30-gauge needle has the smallest diameter.



Pulpdent Code	Description	For Use with	List Price(US\$)
30F50	Orange, 30 ga x 1", straight Tips, pkg. 50		11.40
30F100	Orange, 30 ga x 1", straight Tips, pkg. 100	File-Rite	22.80
25B20	Light Blue, 25 ga x 1/2", Prebent Tips, pkg. 20	Etch-Rite, Etch Royale,	5.00
25B50	Light Blue, 25 ga x 1/2", Prebent Tips, pkg. 50	Etch-All, Porcelain	10.85
25B100	Light Blue, 25 ga x 1/2", Prebent Tips, pkg. 100	Etch Gel	21.75
23R20	Red, 23 ga x 1/2", Prebent Tips, pkg. 20	Pit & Fissure Sealant when a	5.00
23R100	Red, 23 ga x 1/2", Prebent Tips, pkg. 100	smaller applicator is preferred	21.75
22D20	Dark Blue, 22 ga x 1/2", Prebent Tips, pkg. 20	Multi-Cal, Forendo Paste,	5.00
22D100	Dark Blue, 22 ga x 1/2", Prebent Tips, pkg. 100	Dentin Conditioning Gel	21.75
22K20	Black, 22 ga x 1/2", Prebent Tips, pkg. 20	Pit & Fissure Sealant when a metal cannula tip is	5.00
22K100	Black, 22 ga x 1/2", Prebent Tips, pkg. 100	preferred	21.75
22DR15	Dark Blue, 22 ga x 1/2", Prebent Red Dropper Tips, pkg. 20	Dry-Rite, Silane, Polyamide Dentin Sealant	7.55
22DR75	Dark Blue, 22 ga x 1/2", Prebent Red Dropper Tips, pkg. 75		33.15
22KK15	Black, 22 ga x 1/2", Prebent Black Dropper Tips, pkg. 15	Pit & Fissure Sealant when a plastic dropper tip	7.55
22KK75	Black, 22 ga x 1/2", Prebent Black Dropper Tips, pkg. 75	is preferred	33.15
20L20	Pink, 20 ga x 1/2", Prebent Tips, pkg. 20	Kool-Dam, Lime Lite	5.00
20L100	Pink, 20 ga x 1/2", Prebent Tips, pkg. 100	Flows-Rite	21.75
19K20	Black, 19 ga x 1/2", Prebent Tips, pkg. 20	Kool-Dam, Embrace	5.00
19K100	Black, 19 ga x 1/2", Prebent Tips, pkg. 100	WetBond Class V	21.75
18G20	Green, 18 ga x 1/2", Prebent Tips, pkg. 20	Kool-Dam when a larger	5.00
18G100	Green, 18 ga x 1/2", Prebent Tips, pkg. 100	applicator tip is preferred	21.75
P2520	Dark Blue, 25 ga x 1/2", All Plastic Flocked Tips, pkg. 20	Embrace Seal-n-Shine,	6.80
P25100	Dark Blue, 25 ga x 1/2", All Plastic Flocked Tips, pkg. 100	Embrace First-Coat	32.00
FD20	Automix Cartridge Tip + Intra-Oral Tip, pkg. 20 each	Spee-Dee Build-Up	28.50
FSB20	Automix tips	Tuff-Temp Plus 50 mL cartridge	21.15
A20	Automix Syringe Tips, pkg. 20	Embrace Resin Cement, Tuff-Temp Plus, ACTIVA	17.30
A50	Automix Syringe Tips, pkg. 50	5 mL syringe	38.50
A20N1	Automix Tips, clear, with bendable 20-gauge metal cannula – pkg of 20	Supplied with ACTIVA BioACTIVE-BASE/LINER and RESTORATIVE	27.35
A50N1	Automix Tips, clear, with bendable 20-gauge metal cannula – pkg of 50	For use with 2.5mL, 5mL and 10mL 1:1 automix syringes	62.35
AD20T	Automix Tips, clear + short intraoral tips (IOT) – pkg of 20	Supplied with ACTIVA BioACTIVE- RESTORATIVE	27.35
AD50T	Automix Tips, clear + short intraoral tips (IOT) – pkg of 50	For use with 2.5, 5 and 10mL 1:1 automix syringes	62.35
AD20R	Automix Tips, clear + long, narrow intraoral tips (IOR) – pkg of 20	For use with 2.5mL, 5mL and 10mL 1:1 automix syringes	27.35
AD50R	Automix Tips, clear + long, narrow intraoral tips (IOR) – pkg of 50		62.35
AS20	Automix Tips, clear, straight, tapered – pkg of 20	For use with 2.5mL, 5mL and 10mL 1:1 automix syringes	16.40
AS50	Automix Tips, clear, straight, tapered – pkg of 50		37.40
DS05	ACTIVA-SPENSER™ - Dispenser for 5mL 1:1 automix syringes	ACTIVA BioACTIVE-RESTORATIVE	92.50
DS24	Dispenser for 25mL 1:1 automix cartridges	Spee-Dee Build-Up	90.00
DS50	Dispenser for 50mL 1:1 automix cartridges	Spee-Dee Build-Up, Tuff-Temp Plus	90.00